

## OWEAR® ISOLATION GOWN

### Re-order Code and description:

- GN001 – AAMI Level 1, Splash Resistant, Regular, Blue, Elastic cuff
- GN002 – AAMI Level 1, Splash Resistant, Regular, White, Elastic cuff
- GN003 – AAMI Level 1, Splash Resistant, Regular, Yellow, Elastic cuff
- GN001P – AAMI Level 2, \*Impervious, Regular, Blue, Sof-cuff
- GN003P – AAMI Level 2, \*Impervious, Regular, Yellow, Sof-cuff
- GN005P – AAMI Level 2, Impervious, Regular, Blue, Elastic cuff
- GN006P – AAMI Level 2, Impervious, Regular, Yellow, Elastic cuff
- GN009P – AAMI Level 2, Impervious, Polybreath, Regular, White, Elastic cuff
- GN036P – AAMI Level 2, \*Impervious, Extra Large, Blue, Sof-cuff
- GN056P – AAMI Level 2, \*Impervious, Extra Large, Yellow, Sof-cuff

\*Impervious when donned with gloves

### Construction:

- |                   |  |
|-------------------|--|
| Material          | <ul style="list-style-type: none"> <li>- <i>GN001, GN002, GN003</i>: SMS Polypropylene</li> <li>- <i>GN001P, GN003P, GN005P, GN006P, GN009P, GN036P, GN056P</i>: Polyethylene + Polypropylene Blend</li> </ul> <p>Elastic cuff: Latex free elastic<br/>Sof-cuff: Cotton</p>  |
| Bonding & Closure | <ul style="list-style-type: none"> <li>- <i>GN001, GN002, GN003</i>: Stitched neckline and cuff, rest is sonically bonded. Tie at neck and waist.</li> <li>- <i>GN001P, GN003P, GN036P, GN056P</i>: Stitched neckline and cuff, rest is sonically bonded. Velcro fastener at the collar and tie at waist.</li> <li>- <i>GN005P, GN006P, GN009P</i>: Sonically bonded throughout all joins. Tie at neck and waist.</li> </ul> |

**Specification:**

Sizes	Regular	Extra Large
Gown Length	125cm	143cm
Gown Width	150cm	160cm
Sleeve Length with Sof-cuff	68cm	68cm
Sleeve Length with Elastic Cuff	65cm	N/A
Shoulder Width	70cm	80cm

Other - Latex free

Shelf Life - 5 years

Inspection Tolerance -  $\pm 5\%$

**Compliance Standard:**

- Class I non sterile
- ARTG number – 245636
- Splash resistant gowns – AAMI Level 1
- Impervious gowns – AAMI Level 2
  - ISO6529 - *Chemotherapy drugs tested- Cyclophosphamide, 5-Fluorouracil, Cisplatin, Dacarbazine(DTIC), Doxorubicin HCl, Etoposide (Toposar), Paclitaxel (Taxol), Thiotepa and Vincristine Sulphate.*

**Recommended storage:**

Store in a clean, dry environment away from direct sunlight.

**Packaging:**

Dispenser: 10 gowns packed in a polybag

Shipper (Selling Unit): 5 Dispenser

**Intended Purpose:**

Owear® Isolation Gowns are intended to be worn by healthcare providers or visitors to isolate them from the patient or to protect them from a contagious agent that has infected the patient.